


POWER TO SAVE A LIFE

HEARTSTART
ONSITE DEFIBRILLATOR

HEARTSTART ONSITE DEFIBRILLATOR



For the Ordinary Person in the Extraordinary Moment.


PHILIPS



Each year sudden cardiac arrest (SCA) strikes approximately a quarter-million people in the U.S. alone. The majority of these people have no warning, since they show no prior symptoms. And, sadly, fewer than 5% survive, often because emergency medical services cannot reach them in time.

When sudden cardiac arrest strikes, the electrical system of the heart short-circuits, most often causing an abnormal rhythm known as ventricular fibrillation. Lacking proper blood flow, the person loses consciousness, stops breathing, and will die unless promptly treated. CPR (cardiopulmonary resuscitation) can help a person in cardiac arrest, but it alone cannot save lives. A "shock" from a defibrillator – defibrillation therapy – is needed to restore the heart's normal pumping rhythm. A victim's best chance of surviving SCA is to receive that shock within 5 minutes of collapse.





The American Heart Association estimates that an additional 40,000 lives could be saved each year with widespread access to defibrillators.

The Movement to Save More Lives

Over the last 20 years, there has been a widespread effort to move defibrillators into communities where they can be accessed and used by trained citizens, who might be present at the onset of SCA or first on the scene. To that end, Philips produces easy-to-use automated external defibrillators. Our HeartStart defibrillators enable almost anybody to treat SCA quickly and effectively, wherever it happens — at work, at play, in the air — providing the power to save a life.

The HeartStart OnSite Defibrillator

Philips, the leader in defibrillation technology,* designed the HeartStart OnSite Defibrillator for the ordinary person in the extraordinary moment. Our innovative technology, based on extensive research and user feedback, has produced a defibrillator so easy to use that, with minimal training, you can potentially save the life of a co-worker, friend, or anyone else stricken with sudden cardiac arrest.

Weighing just 3.3 lbs., this small and lightweight defibrillator can be easily carried to the patient's side. Using clear, natural voice instructions, the OnSite Defibrillator guides you, the first caregiver on the scene, through each step of defibrillation and CPR. Integrated SMART Pads placed on the patient's bare skin transmit information to the defibrillator, which senses and adapts to your actions every step of the way.

HeartStart OnSite includes proven Philips technologies for heart rhythm assessment (SMART Analysis) and defibrillation energy delivery (SMART Biphasic). And like all HeartStart defibrillators, it can be used to treat infants, children and adults.

Ready When Needed

HeartStart Onsite is powered by a simple, safe, long-life (four-year) battery, the same technology used in cameras. Just click the battery in place, and the defibrillator is standing by, ready for use.

Between uses, the defibrillator performs daily, weekly, and monthly self-tests of its electrical components, subsystems, and battery to help ensure continued reliability and readiness. A pads integrity test, which checks daily that the cartridge is properly installed and that the defibrillator pads are in working order, provides an additional level of assurance.

As long as the green "Ready" light on the defibrillator is blinking, you can be certain the HeartStart OnSite is ready for use. If any part of the defibrillator needs attention, the "Ready" light turns off and the defibrillator "chirps" to alert you. Pressing the flashing blue "Information" button provides verbal guidance to help you identify and correct the situation.

*Frost & Sullivan

HeartStart OnSite Defibrillator is designed to enable you to help save a life in extraordinary circumstances.

- ▶ **Lightweight:** Fully equipped at just 3.3 pounds.
- ▶ **Intuitive:** Clean design and clear voice instructions instill the confidence that's needed when administering therapy to a person in cardiac arrest.
- ▶ **Ever-Ready:** Powered by a long-life (4-year) disposable battery. Automated tests help ensure readiness.
- ▶ **Versatile:** Available for use on anyone of any age, including children and infants.
- ▶ **Effective:** Patented SMART Analysis heart rhythm assessment and SMART Biphasic defibrillation therapy, clinically proven in nearly 10 years of use. No other biphasic waveform is as well documented.



Easy To Use

Using the HeartStart OnSite Defibrillator is simple. Pulling the green handle powers-up the defibrillator and activates voice instructions. These instructions are paced to your actions, to help guide you through the entire process, from placing pads on the patient to delivering a defibrillation shock.


HeartStart OnSite determines if a shock is needed, based on its automated assessment of the victim's heart rhythm:

- ▶ If a shock is needed, the defibrillator directs you to press the flashing orange "Shock" button. Then, HeartStart OnSite delivers a dose of low-energy biphasic therapy, a highly effective defibrillation waveform that is also gentle to the heart.
- ▶ If a shock is not needed, the OnSite Defibrillator instructs you to assess the patient and to perform CPR when necessary. While performing CPR, the defibrillator's voice instructions can be activated to coach you on the frequency and depth of compressions.

HeartStart OnSite also reminds you to call emergency medical services (EMS). And, should EMS need a summary of care, it can be retrieved from the defibrillator's internal memory. An EMS provider simply presses the "i-button" and HeartStart Onsite verbally recounts events from its last clinical use.

Replaceable SMART Pads Cartridges

The cartridge contains two adhesive pads that are peeled and placed on the patient's bare skin as indicated by the pictures on the pads. The pads are "smart" because they can sense that they have been removed from the cartridge and applied to the patient.



HeartStart OnSite can be used on patients of any age. Special Infant/Child SMART Pads reduce the defibrillator's shock energy, so that you can treat an infant or child suffering sudden cardiac arrest. The Infant/Child Pads cartridge is marked with an indication of the appropriate weight and with a teddy bear icon for ready identification.

To practice your skills, a special training pads cartridge (adult or infant/child) can be installed in the defibrillator. It suspends the defibrillator's ability to shock, while walking you through patient care scenarios.

DEFIBRILLATOR

Defibrillator Model	HeartStart M5066A
Defibrillator Family	HS1
How Supplied	Defibrillator, Instructions for Use, battery, adult SMART Pads cartridge (2), Quick Reference
Waveform	Truncated Exponential Biphasic. Waveform parameters adjusted as a function of each patient's impedance.
Energy	Single energy output. Adult: nominal 150 Joules into a 50 ohm load. Infant/Child: nominal 50 Joules into a 50 ohm load. Automatically set based on type of SMART Pads cartridge installed.
Shock-to-Shock Cycle Time	Typically less than 20 seconds.
Protocol	Voice prompts and indicators guide user through protocol. Follows preconfigured settings. Can be modified with HeartStart Event Review software.
Voice Instructions	Detailed voice messages guide responder through use of the defibrillator.
CPR Voice Coaching	Voice instructions in adult and infant/child CPR available at user's option.
Shock Delivery	Via adhesive pads placed on patient's bare skin as illustrated on pads.
Controls	Green SMART Pads cartridge handle, Green On/Off button, Blue i-button, Orange shock button.
Indicators	Ready light; Blue i-button; Caution light.

PHYSICAL

Size	
Height	2.8" (7 cm)
Width	8.3" (21 cm)
Depth	7.4" (19 cm)
Weight	
With battery and pads cartridge	3.3 lbs. (1.5 kg)
Without battery or pads cartridge	2.4 lbs. (1 kg)

ENVIRONMENTAL/PHYSICAL REQUIREMENTS

Sealing	Drip proof per EN60529 class IPX1. Solid Objects per EN60529 class IP2X.
Temperature	Operating: 32° - 122° F (0° - 50° C) Standby: 50° - 109° F (10° - 43° C)
Humidity	Operating: 0% to 95% relative, non-condensing Standby: 0% to 75% relative, non-condensing
Altitude	Operating: 0 to 15,000 feet Standby: 0 to 8,500 feet >48 hours and 8,500 to 15,000 feet <48 hours
Shock/Drop Abuse	Withstands 1 meter drop to any edge, corner or surface.
Vibration	Meets EN1785 random and swept sine, road ambulance specification in operating and standby states.
EMI (Radiated/Immunity)	Meets EN55011 Group 1 Level B Class B and EN61000-4-3.

PATIENT ANALYSIS SYSTEM

Patient Analysis	Evaluates patient ECG and signal quality to determine if a shock is appropriate, and evaluates connection impedance for proper defibrillation pad contact.
Sensitivity/Specificity	Meets AAMI DF39 guidelines and AHA recommendations for adult defibrillation (<i>Circulation</i> 1997;95:1677-1682).
Artifact Detection	The effects of pacemaker artifact and electrical noise are minimized with artifact detection.

BATTERY (M5070A)

Type	9 Volt DC, 4.2 Ah, lithium manganese dioxide, disposable long-life primary cell.
Capacity	Minimum 90 shocks or 3 hours of operating time.
Install-By Date	Battery is labeled with an install-by date of at least five years from date of manufacture.
Standby Life	Four years typical when battery is installed by the install-by date. (Will power the AED in standby state within the specified standby temperature range, assuming one battery insertion test and no defibrillation uses.)

SMART PADS

Adult SMART Pads Cartridge	M5071A defibrillation pads for patients 8 years of age and older or 55 lbs. (25 kg) and over.
Infant/Child SMART Pads Cartridge	M5072A defibrillation pads for patients under 8 years of age or 55 lbs. (25 kg).
Energy Delivered	Adult: nominal 150 Joules into a 50 ohm load. Infant/Child: nominal 50 Joules into a 50 ohm load.
How Supplied	Disposable cartridge, containing adhesive defibrillation pads, clicks into defibrillator for an integrated pads solution.
Active Surface Area	13.2 in ² (85 cm ²) each
Cable Length	Adult pads: 54 inch (137.1 cm) Infant/Child pads: 40 inch (101.6 cm)
Use-By Date	Cartridge is labeled with a use-by date of at least two years from date of manufacture.

TRAINING PADS

Adult Training Pads Cartridge	M5073A
Infant/Child Training Pads Cartridge	M5074A
Function	Special pads place HeartStart OnSite into training mode and disable its energy delivery capability. Feature 8 real-world training scripts.

AUTOMATED AND USER-ACTIVATED SELF-TESTS

Daily Automatic Self-tests	Tests internal circuitry, waveform delivery system, pads cartridge and battery capacity.
Pads Integrity Test	Specifically tests readiness-for-use of pads (gel moisture).
Battery Insertion Test	Upon battery insertion, extensive automatic self-tests and user-interactive test check device readiness.
Status Indicator	Blinking green Ready light indicates ready for use.

DATA RECORDING AND TRANSMISSION

Infrared	Wireless transmission of event data to personal computer using the IrDA protocol.
HeartStart Event Review Software	Data management software "optional" for download and review of data retrieved through defibrillator's infrared data port.
Data Stored	First 15 minutes of ECG and the entire incident's events and analysis decisions.

Refer to HeartStart OnSite Defibrillator Instructions for Use for detailed product instructions. All specifications based on 25° C unless otherwise noted. The defibrillator and its accessories are made of latex-free materials.

Philips, Leading the Way

Philips, the leader in defibrillation technology, offers a portfolio of resuscitation products to help rescuers treat SCA quickly and effectively wherever it occurs.

Philips automated defibrillators are located in a wide range of settings, such as office buildings, industrial complexes, hospitals, airports and airplanes, resorts, fitness centers, shopping malls and educational facilities.

As the newest addition to the Philips line of defibrillators, HeartStart OnSite optimizes advanced defibrillation technology and makes it available to a broad range of users. Now, this therapy can be seconds away in the home, office, or in any public or private location where people work, live and play.

To learn more about HeartStart OnSite Defibrillator and Philips Medical Systems, visit www.medical.philips.com/cms or call 800-453-6860.



PHILIPS MEDICAL SYSTEMS IS PART OF ROYAL PHILIPS ELECTRONICS

On the web: www.medical.philips.com

Via e-mail: medical@philips.com

By fax: +31 40 27 64 887

By postal service:

Philips Medical Systems
 Global Information Center
 I.B.R.S. / C.C.R.I. Numéro 11088
 5600 VC Eindhoven
 Pays-Bas / The Netherlands

CORPORATE HEADQUARTERS

Philips Medical Systems
 Nederland B.V.
 Postbus 10.000
 5680 DA Best
 Netherlands

UNITED STATES

Andover, Massachusetts
 (800) 934-7372

CANADA

Markham, Ontario
 (800) 291-6743

EUROPE, MIDDLE EAST AND AFRICA

Böblingen, Germany
 Fax: (49) 7031 463 1552

LATIN AMERICA HEADQUARTERS

Sunrise, Florida
 (954) 835-2600

ASIA PACIFIC HEADQUARTERS

Wanchai, Hong Kong
 (852) 2821 5888

© Philips Medical Systems
 North America Corporation, 2002
 All rights are reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright holder.

Philips Medical Systems North America Corporation reserves the right to make changes in specifications or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Printed in U.S. November 2002
 Part #5990-0572EN

NOTE: The HeartStart OnSite Defibrillator and some of its accessories are available with a prescription.